EPI-NO *Libra*

Pelvic floor trainer

Beckenbodentrainer

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Product Description

The **EPI-NO Libra** Pelvic Floor Muscle Trainer was designed, developed and manufactured in Germany in close co-operation with gynaecologists, physiotherapists and women.

It consists of (1) a narrow contoured silicone balloon, (2) a hand pump, (3) a pressure display, (4) an air release valve, (5) connected by a flexible plastic tube.

EPI-NO Libra

Indication

EPI-NO Libra is patented and intended to be used as a Pelvic Floor Muscle Trainer.

Application

EPI-NO Libra promotes optimal pelvic floor muscle training.

Contraindications

You must not use **EPI-NO Libra** if any of the following are present:

- There is a risk of an ascending infection. (e.g. vaginal fungal infection).
- You have injuries in your genital area which have not yet healed.
- You are under the influence of alcohol or drugs.
- You have a malignant disease in your genital area (e.g. cervical cancer)

EPI-NO Libra should only be used under the supervision of a physician or physiotherapist:

- If you have any nerve damage or other diseases that reduce pain sensation in your genital area (e.g. permanently high blood sugar levels).
- If you are a paraplegic
- If you have been diagnosed with Multiple Sclerosis.
- If pain killers that reduce the feeling of pain in the vaginal area have been administered.
- If you have marked varicosity at the vaginal entrance (varicose veins).
- If you have adhesions or abnormalities in the vaginal or cervicall area

Please consult your physician or physiotherapist before using **EPI-NO Libra**.

Warning / Caution

- For hygienic reasons, **EPI-NO Libra** is intended as a single user device; meaning only one woman should use the device.
- Do not use **EPI-NO Libra** after the end of its recommended life span (as shown on the package label of the Balloon Tube System). It is recommended, for hygiene reasons, to replace the balloon-tube-system 6 months after first use.
- Stop exercising with **EPI-NO Libra** if you experience any pain.
- Before you use **EPI-NO Libra** the first time or if it has not been used for more than four weeks, you should inflate the balloon by first closing the air release valve, then squeezing the hand pump about ten times in sequence. Open the air release valve to let the air out and return the balloon to its original shape. This allows you to test the device. If your **EPI-NO Libra** is not in perfect working order contact your supplier before use.
- We recommend you thoroughly wash your hands with soap and water before each use of **EPI-NO Libra**. Please strictly observe the guidelines for cleaning and disinfecting **EPI-NO Libra** (see section "Cleaning and Disinfecting").
- Only use water-based lubricants. Do not use lubricants containing oils or essential components as these can damage the balloon.
- Take care not to damage the balloon with pointy or sharp-edged items (e.g. rings, fingernails, etc)
- **EPI-NO Libra** must not be used under water (e.g. in the bath) as penetrating water can impair function.
- Store **EPI-NO Libra** carefully to avoid extended exposure to direct sunlight.
- **EPI-NO Libra** should be kept out of reach of children!


The Function of the Pelvic Floor

The pelvic floor consists of a sling of muscles stretching from the tailbone at the base of your spine to the pubic bone at the front. This important muscle group supports the uterus, bladder and bowel.

A strong pelvic floor supports the stability of the spinal column and positively affects the entire body posture. With its mesh-like connections, it seals the bottom of the pelvis.

It is easy to understand that weakness in this muscle group can lead to sagging of the bladder, uterus and vagina leading to unpleasant side effects:

These can range from impaired sexual function to incontinence.

Pelvic floor muscles (red)

Typical symptoms of weakened pelvic floor muscles are:

- Leaking of urine when coughing, sneezing, climbing stairs, and carrying heavy objects.
- Constant urge to urinate, even with an empty bladder.

- A strong feeling of downward pressure.
- Prolapse of the vaginal tissues or the uterus.
- Involuntary loss of wind (flatus) and stool (anal incontinence).

Pelvic Floor Muscle Exercise

Select a pleasant place for your exercise with **EPI-NO Libra**. Ensure that you will not be disturbed.

Choose a comfortable position (e.g. semi reclined). With time and practice you will find the most appropriate position for you.

Take relaxed, regular breaths during exercise. Try to relax your buttocks, thighs and stomach. We recommend exercising 10 to 20 minutes per day.

Getting Started

Open the air release valve, squeeze the air out of the balloon with your hand, then close the air release valve.

The balloon will now remain in this deflated shape. Lubricating the top half of the balloon will facilitate insertion (see section "Warning / Caution").

Hold the balloon near the end and introduce it about two thirds of the way into the vagina. Give yourself time: The vagina may be dry and sensitive.


Relax your pelvic floor. Now inflate the balloon to the point that you are aware of its presence (for most women the pressure gauge will show 2).

Should the pressure be too strong, you can reduce it at any time by opening the air release valve.

Perineal Draw Up Exercise

Keep a diary of your activity.

Exercise sessions should occur twice a day until you have reached your treatment goal, there after once per day for maintenance.

At the start of your session you should note your resting vaginal pressure on the **EPI-NO Libra** pressure gauge.

While tensing the pelvic floor muscles try to draw the balloon upward for about 10 seconds (in the beginning your pelvic floor muscles may be weak and you can build up to 10 seconds). Then relax for 10 seconds and repeat this exercise again.

Gently touch the balloon and feel it moving upwards with the muscle tension. This way you can check that you are exercising correctly.

Observe the needle on the pressure display before and after tensing your pelvic floor muscles.

The difference between resting pressure and contraction pressure indicates the contraction strength: See "Training Progress Table".

You may not see any appreciable pressure difference when you first commence your exercises. However, after a few days you will begin to achieve your first results which will show on the **EPI-NO Libra** pressure gauge.

Ideally, pelvic floor muscle exercises should be continued life-long (see section "Accessor/consumable article").

If any complications should arise when using the **EPI-NO Libra** please consult your clinician immediately.

The further exercises below may be added for variety.

Further exercises with EPI-NO Libra

The following exercises give you a variety of training possibilities. Follow the instructions for use or train with the **EPI-NO Libra** exercises recommended by your clinician.

Follow the instructions for inserting your **EPI-NO Libra** (see section "Pelvic Floor Exercises").

During all exercises you can follow the needle on the pressure gauge and check your pelvic floor strength.

Training Progress Table:

Pressure difference (tensed - relaxed)	Comment
< 1	Should be improved
1 - 2	You are getting there
3 - 4	Already very good
> 4	Excellent

“Pelvis Rock”

Lie on your back and bend your knees with your feet firmly on the ground, and pressure on your heels. Alternately raise your right and left hip and repeat this 2-3 times.

“Around the Clock”

Assume the same starting position as for “Pelvis Rock”. Imagine there is a clock under your pelvis.

Follow the movement of the hands, first clockwise and then counter clockwise. The movements should be so smooth that the pelvis seems to continuously rotate.

“Long Legs”

Lie on your back and stretch your legs. The entire spinal column should touch the floor. Move your chin towards your chest so that the neck is stretched.

Alternately push your right and left heel away as if you want to elongate your legs. Simultaneously draw your toes towards your body.

Cleaning and Disinfecting

Caution

- Make sure the pump does not get wet during cleaning and disinfecting because water may damage the interior.
- Follow cleaning and disinfecting instructions prior to first use and immediately after every subsequent use.
- If **EPI-NO Libra** is not used for more than 7 days, repeat cleaning and disinfecting prior to use.

Before first usage and immediately after every use, clean **EPI-NO Libra** thoroughly following these guidelines :

- Inflate the balloon by pressing the pump about 10 times.
- Carefully clean the balloon with soap and warm water (minimum 40°C / 104°F; 1-2 minutes). Take care not to wet the pump as water may damage the interior.
Rinse the balloon under warm running water. Let **EPI-NO Libra** dry on a clean, smooth surface.
- Disinfect **EPI-NO Libra** by wiping it with 70% Isopropyl from your pharmacy. Use a cotton pad or medical alcohol pads.
- Only when completely dry put **EPI-NO Libra** into the pouch, close zip and store in a clean, dry place.

We wish you success
with your **EPI-NO Libra** training.

Accessor / consumable article

It is recommended, for hygienic reasons, to replace the Balloon-Tube-System (**EPI-NO Libra** BTS) every 6 months, if used regularly.

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TECSANA GmbH